

BioXp™ 3200 system

The BioXp™ 3200 system is an automated workstation for a variety of genomic and synthetic biology workflows, including generating clones, DNA variant libraries, genes, and genomes.

Revolutionizing custom DNA synthesis

Say hello to the BioXp™ system and say goodbye to waiting weeks for service labs or synthesizing genes on the bench.

The BioXp™ system can build up to 32 complex constructs, from 300 to 7,000 base pairs, with multiple fragments and vectors, in less than 21 hours. Bypass plasmid preps and synthesize transfection-ready DNA quantities of more than 8 µg. The BioXp™ system allows you to eliminate cloning and synthesize *de novo* gene constructs with sequence modifications, accelerating the development of vaccines and biologics.



Build better. Build faster. Build cheaper. Only on the BioXp™ system.



Control your timelines

Free yourself from delayed orders – BioXp™ application orders arrive within 3 to 5 days



Maximize your lab budget

BioXp™ applications are value-priced and come with everything you need



Secure your intellectual property

Enjoy peace of mind, knowing that your vectors never leave your lab

BioXp™ applications

Building biology on your lab bench doesn't get easier than with the BioXp™ system. From building gene fragments and clones to library construction, BioXp™ applications come with everything you need to complete projects.



Build genes

Custom high-fidelity gene fragments



Build clones

Custom DNA fragments, up to 7,000 bp, cloned into any vector




Build libraries

Custom variant DNA libraries

Design. Build. Discover.

Enjoy the freedom of controlling your timelines with the BioXp™ system. Place BioXp™ application orders using our online portal, and everything you need is shipped to you within three to five days. Load the instrument with the provided labware and press *start*. Collect your constructs the next day. It's as simple as that.



Building biology has never been easier.



Order DNA

Log in to the BioXp™ ordering portal and submit your DNA sequences



Receive application kit

Kits containing the pieces of your desired gene(s) arrive in three to five days*



Push *start* and walk away

Load the cartridges and walk away — your finished product can be ready in as little as six hours

Product specifications

Weight	63.4 kg [139.8 US lbs]
Dimensions (W × D × H)	66 × 74 × 53 cm [26 × 29 × 21 inches]

*International delivery times may vary.
BioXp™ and the Codex DNA logo are registered trademarks or trademarks of Codex DNA, Inc.
For research use only. Not for use in diagnostic procedures.

Codex DNA, Inc.
9535 Waples Street, Suite 100
San Diego, CA 92121 USA
hello@codexdna.com | **858.228.4115** | codexdna.com